



Recombinant Nautilia profundicola Protein CrcB homolog (crcB)

Product Code	CSB-CF497834NFR
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B9L8U9
Product Type	Transmembrane Protein
Immunogen Species	Nautilia profundicola (strain ATCC BAA-1463 / DSM 18972 / AmH)
Sequence	MKFDTMLAIGIGGFIGAILRAYTAGLVNSAVKHDIPFGTLSVNLIGSLLLGMFIGAI QYG GIQNPYIKSMLTTGMMGAFTTFSTFAVESFFLTKNALYIQALS YILLNVIGCIILAG AGF KGIEAILK
Source	in vitro E.coli expression system
Target Names	crcB
Protein Names	Recommended name: Protein CrcB homolog
Expression Region	1-128
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.