



Recombinant Arabidopsis thaliana Ergosterol biosynthetic protein 28 (At1g10030)

Product Code	CSB-CF526812DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O80594
Product Type	Transmembrane Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Sequence	MKALGYWLMVVGSLRLASVWFGFFNIWALRLAVFSQTTMSEVHGRTFGVWTL LTCTLCFL CAFNLNKPPLYLATFLSFIYALGHFLTEYLFYQTM TIANLSTVGFFAGTSIVWML LEWNS LEQPHSKLS
Source	in vitro E.coli expression system
Target Names	At1g10030
Protein Names	Recommended name: Ergosterol biosynthetic protein 28
Expression Region	1-129
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.