



Recombinant Staphylococcus carnosus UPF0344 protein Sca_0577 (Sca_0577)

Product Code	CSB-CF494229FLK
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B9DIS4
Product Type	Transmembrane Protein
Immunogen Species	Staphylococcus carnosus (strain TM300)
Sequence	MLHLHIFSWVIGIILFIVSYISFTKTGAPKKAYKPLHMTLRLFLVLILFSGVWQVVE EFA TATGSTHMLLTLKMICGIGVVALMEVTLVRKQRGASHKGLFWGTIALIIVTMALG IILPG GPISNMFVITK
Source	in vitro E.coli expression system
Target Names	Sca_0577
Protein Names	Recommended name: UPF0344 protein Sca_0577
Expression Region	1-131
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.