



Recombinant Human Leptin receptor gene-related protein (LEPROT)

Product Code	CSB-CF012875HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O15243
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MAGVKALVALSFSGAIGLTFLMLGCALEDYGVYWPLFVLIFHAISPIPHFIKRVT YDSD ATSSACRELAYFFTTGIVVSAFGFPVILARVAVIKWGACGLVLAGNAVIFLTIQG FFLIF GRGDDFSWEQW
Source	in vitro E.coli expression system
Target Names	LEPROT
Protein Names	Recommended name: Leptin receptor gene-related protein Alternative name(s): Leptin receptor overlapping transcript protein OB-R gene-related protein Short name= OB-RGRP
Expression Region	1-131
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.