



Recombinant Escherichia coli O81 UPF0299 membrane protein yohJ (yohJ)

Product Code	CSB-CF485342EOP
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B7MXF3
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli O81 (strain ED1a)
Sequence	MSKTLNIIWQYLRAFVLIYACLYAGIFIASLLPVTIPGSIIGMLILFVLLALQILPAKW V NPGCYVLIRYMALLFVPIGVGVMQYFDLLRAQFGPVVVSCAVSTLVVFLVVS W SSQLVHG ERKVVGQKGSEE
Source	in vitro E.coli expression system
Target Names	yohJ
Protein Names	Recommended name: UPF0299 membrane protein yohJ
Expression Region	1-132
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.