



Recombinant Escherichia coli O139:H28 Protein psiE (psiE)

Product Code	CSB-CF422004EJD
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A7ZUQ1
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli O139:H28 (strain E24377A / ETEC)
Sequence	MTSLSRPRVEFISTILQTVLNLGLLCLGLILVVFLGKETVHLADVLFAP EQTSKYE LVEG LVVYFLYFEFIALIVKYFQSGFHFPLRYFVYIGITAIVRLIIVDHKSP LDVLIYSAAIL L LVITLWLCNSKRLKRE
Source	in vitro E.coli expression system
Target Names	psiE
Protein Names	Recommended name: Protein psiE
Expression Region	1-136
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.