



Recombinant Human Transmembrane protein 203 (TMEM203)

Product Code	CSB-CF023797HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q969S6
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MLFSLRELVQWLGFATFEIFVHLLALLVFSVLLALRVDGLVPGLSWWNVFVPFF AADGLS TYFTTIVSVRLFQDGEKRLAVLRLFWVLTVLSLKVFVFEMLLCQKLAEQTRELWF GLITSP LFILLQLLMIRACRVN
Source	in vitro E.coli expression system
Target Names	TMEM203
Protein Names	Recommended name: Transmembrane protein 203
Expression Region	1-136
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.