



Recombinant Pasteurella multocida UPF0299 membrane protein PM0880 (PM0880)

Product Code	CSB-CF875062ESG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9CME9
Product Type	Transmembrane Protein
Immunogen Species	Pasteurella multocida (strain Pm70)
Sequence	MTRKIVDLARSCGILYLMLFIGEWIAHHLNIGIPASIWGLLLLFLGLTFRIIKLDWV LCS ASLLIRYMALLFVPVSVGVIKYADVLFSQMNVLLLPNIVSTFLTLIVVGLLSDYLF SLSS FSHLRKKVARKQAEKA
Source	in vitro E.coli expression system
Target Names	PM0880
Protein Names	Recommended name: UPF0299 membrane protein PM0880
Expression Region	1-136
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.