



Recombinant Inner membrane protein yphA (yphA)

Product Code	CSB-CF359885SZB
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0AD48
Product Type	Transmembrane Protein
Immunogen Species	Shigella flexneri
Sequence	MNTRLRYFDFGAARPVLLLIARIAVVLIFIFGFPMKMMGFDGTVQYMASLGAPMP MLAAII AVVMEVPAAILIVLGFFTRPLAVLFIFYTLGTAVIGHHYWDMTGDAVGPNMINF WKNVSI AGAFLLLAITGPGAISLDRR
Source	in vitro E.coli expression system
Target Names	yphA
Protein Names	Recommended name: Inner membrane protein yphA
Expression Region	1-140
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.