



Recombinant Bacillus subtilis Uncharacterized protein yxaJ (yxaJ)

Product Code	CSB-CF331429BRJ
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P42109
Product Type	Transmembrane Protein
Immunogen Species	Bacillus subtilis (strain 168)
Sequence	MSKKNSVAIMTTISAFLEFCVIVAASLSPLAGTGAAANQFNSAGMWSAVGMILV LYFIPF LVYMLGVDAMRYVMAVLCGFGLLIHLSSAGFILMFSFFSDHLLSEVIFVLGVSLA AAAVN VIWFVAAFVRSAGAEKTSVNTLT
Source	in vitro E.coli expression system
Target Names	yxaJ
Protein Names	Recommended name: Uncharacterized protein yxaJ
Expression Region	1-141
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.