



# Recombinant *Ashbya gossypii* Golgi apparatus membrane protein TVP15 (TVP15)

<b>Product Code</b>	CSB-CF748324DOT
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	Q74ZU2
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	<i>Ashbya gossypii</i> (strain ATCC 10895 / CBS 109.51 / FGSC 9923 / NRRL Y-1056) (Yeast) ( <i>Eremothecium gossypii</i> )
<b>Sequence</b>	MSENDQKGVKFIRALLVTTGAVATLGLVQLSEVASEFLAAARACFGLPLSVLL VYLEFR PVPLLQQYASFYYSYMGRALLQLLLAVMLLPAGFFQVCAFTMLFMTGFICAALE LSSSAP ELPSFRNDGRALSIGAEEDDDMI
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	TVP15
<b>Protein Names</b>	Recommended name: Golgi apparatus membrane protein TVP15
<b>Expression Region</b>	1-143
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.