



Recombinant *Listeria innocua* serovar 6a Cell wall teichoic acid glycosylation protein gtcA (gtcA)

Product Code	CSB-CF842158LGU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q927U4
Product Type	Transmembrane Protein
Immunogen Species	<i>Listeria innocua</i> serovar 6a (strain ATCC BAA-680 / CLIP 11262)
Sequence	MSRVRQLLDKLPWYTDNIHSILMYLIMGGFTTLINIVTFWLCTDVLNWDYRIANT IAWVA SVLFAYFSNKKYVFESYPTWKEKAREVSSFFGFRFLTYIVDFLVMILLISGLGIN ELWA KIWTNVIVLILNRYVFSKWIIFKVRK
Source	in vitro E.coli expression system
Target Names	gtcA
Protein Names	Recommended name: Cell wall teichoic acid glycosylation protein gtcA
Expression Region	1-145
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.