



Recombinant Uncharacterized protein Rv2081c/MT2143 (Rv2081c, MT2143)

Product Code	CSB-CF611995MVZ
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q10689
Product Type	Transmembrane Protein
Sequence	MFANAGLSPFVAIWTARAASLYTSHNFWCAA AVSAAVYVGS AVVPA AVAGPL FVGRVSAT IKAAAPSTTAAIATLATAANGQLRERGGAGGWVGVHCPVVG GGGVGHPRKAI AAAVSVHS TCMPAAF GGHLGLGDRSR SVLSGTP
Source	in vitro E.coli expression system
Target Names	Rv2081c, MT2143
Protein Names	Recommended name: Uncharacterized protein Rv2081c/MT2143
Expression Region	1-146
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.