



Recombinant Escherichia coli UPF0756 membrane protein YeaL (yeaL)

Product Code	CSB-CF471824ENW
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B6IBL2
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli (strain SE11)
Sequence	MFDVTLILLGLAALGFISHNTTVAVSILVLIIVRVTPPLSTFFPWIEKQGLSIGIIILTIGVMAPIASGTLPPSTLIHSFLNWKSLVAIavgVIVSWLGGRGVTLMGSQPQLVAGLLVGT VLGVALFRGVPVGPLIAAGLVSLIVGKQ
Source	in vitro E.coli expression system
Target Names	yeaL
Protein Names	Recommended name: UPF0756 membrane protein YeaL
Expression Region	1-148
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.