



Recombinant Escherichia coli O157:H7 UPF0756 membrane protein YeaL (yeaL)

Product Code	CSB-CF509760EOF
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	C6UX49
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli O157:H7 (strain TW14359 / EHEC)
Sequence	MFDVTLILLGLAALGFISHNTTVAVSILVLIIVRVTPPLSTFFPWIEKQGLSIGIIILTI GVMAPIASGTLPPSTLIHSFLNWKSLVAIavgVIVSWLGGRGVTLMGSQPQLVA GLLVGT VLGVALFRGVPVGPLIATGLVSLIVGKQ
Source	in vitro E.coli expression system
Target Names	yeaL
Protein Names	Recommended name: UPF0756 membrane protein YeaL
Expression Region	1-148
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.