



Recombinant Escherichia coli O6:K15:H31 UPF0756 membrane protein YeaL (yeaL)

Product Code	CSB-CF610681EGY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q0TH40
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
Sequence	MFDVTLILLGLAALGFISHNTTVAVSILVLIIVRVTPPLSTFFPWIEKQGLSIGIIILTI GVMAPIASGTLPPSTLIHSFLNWKSLVAIavgVIVSWLGGRGVTLMGSQPQLVA GLLVGT VLGVALFRGVPVGPLIAAGLVSLIVGKQ
Source	in vitro E.coli expression system
Target Names	yeaL
Protein Names	Recommended name: UPF0756 membrane protein YeaL
Expression Region	1-148
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.