



Recombinant Arabidopsis thaliana CASP-like protein At1g49405 (At1g49405)

Product Code	CSB-CF660368DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q3ECT8
Product Type	Transmembrane Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Sequence	MVEVPGSVGTTASLSLRLGQMVLAFGSLLFMTIGVRFYQFTAFCYLVTIMSLAI PWNLT AMVDIYCVILQQPFQKPRILLAISIGDWVVSVLALASASSAASVVDILRSNESSC PPTIC NRYQFAATLAFITWFLSLSSSLFNLWLLPSLI
Source	in vitro E.coli expression system
Target Names	At1g49405
Protein Names	Recommended name: CASP-like protein At1g49405
Expression Region	1-152
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.