



Recombinant Arabidopsis thaliana CASP-like protein At3g23200 (At3g23200)

Product Code	CSB-CF810541DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8L7R5
Product Type	Transmembrane Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Sequence	MIDIPGTPGTLTGLVLRISQCVFAAGSISYMTSGGFFSFTAFCYLIAAMGLQVI WSFGL AILDTFALVRKKTLLSPVLVSLFVVGDWVTSTLSLAGASSSAGITVLYFGDLGSC SFEAE CWKYQLSVALAFLCWITIAVSSLTTLWLLASA
Source	in vitro E.coli expression system
Target Names	At3g23200
Protein Names	Recommended name: CASP-like protein At3g23200
Expression Region	1-152
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.