



Recombinant *Macaca fascicularis* NADH-ubiquinone oxidoreductase chain 4 (MT-ND4)

Product Code	CSB-CF015079MOV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P50664
Product Type	Transmembrane Protein
Immunogen Species	<i>Macaca fascicularis</i> (Crab-eating macaque) (<i>Cynomolgus</i> monkey)
Sequence	SFSGATTLMIAHGLTSSMYFCLANSNYERTHNRTMLLSRGLQILLPLTAFWWLT ASLTNL ALPPTINLLGELFVITTSFSWSHITIVLTGLNMLITALYSLHMFITVQRGTLTHHMI NMK PPF TRENMLMFMHLAPIILLSLNPNIILGFTS
Source	in vitro <i>E.coli</i> expression system
Target Names	MT-ND4
Protein Names	Recommended name: NADH-ubiquinone oxidoreductase chain 4 EC= 1.6.5.3 Alternative name(s): NADH dehydrogenase subunit 4
Expression Region	1-152
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.