



Recombinant Human Myelin and lymphocyte protein (MAL)

Product Code	CSB-CF013372HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P21145
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MAPAAATGGSTLPSGFSVFTLLPDLLFIFEFFIFGGLVWILVASSLVPWPLVQGW VMFVSV FCFVATTTLLIILYIIGAHHGGETSWVTLDAAAYHCTAALFYLSASVLEALATITMQDG FTYR HYHENIAAVVFSYIATLLYVVHAVFSLIRWKSS
Source	in vitro E.coli expression system
Target Names	MAL
Protein Names	Recommended name: Myelin and lymphocyte protein Alternative name(s): T-lymphocyte maturation-associated protein
Expression Region	1-153
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.