



Recombinant *Macaca mulatta* (Rhesus macaque) Arachidonate 5-lipoxygenase-activating protein (ALOX5AP)

Product Code	CSB-CF001625MOW
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P30354
Product Type	Transmembrane Protein
Immunogen Species	<i>Macaca mulatta</i> (Rhesus macaque)
Sequence	MDQETVGNVLLAIIVTLISVVQNGFFAHKVEHESRTQNGRSFQRTGTLAFERV YTANQNC VDAYPTFLAVLWSAGLLCSQVPAAFAGLMYLLVRQKYFVG YLGERTQSTPGYI FGKRIIL FLFLMSVAGIFNYLIFFFGSDFENYIKTVTTT
Source	in vitro <i>E.coli</i> expression system
Target Names	ALOX5AP
Protein Names	Recommended name: Arachidonate 5-lipoxygenase-activating protein Alternative name(s): FLAP MK-886-binding protein
Expression Region	1-153
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.