



Recombinant Exiguobacterium sp. UPF0756 membrane protein EAT1b_0668 (EAT1b_0668)

Product Code	CSB-CF504714EPU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	C4L477
Product Type	Transmembrane Protein
Immunogen Species	Exiguobacterium sp. (strain ATCC BAA-1283 / AT1b)
Sequence	MSAYLFLMLVLIGVISNNQSVIIASSVLLIIKAIGFGDQLFPTLASKGISWGVTIITI A VLVPIATGDIGFKELWNSIKGPVGIVAFASGMFVAIAAGQGVQLMRVDPVVTTA LLAGTI LAVGFMKGIPVGPLVGAGIAALILGGYQVIEKWF
Source	in vitro E.coli expression system
Target Names	EAT1b_0668
Protein Names	Recommended name: UPF0756 membrane protein EAT1b_0668
Expression Region	1-154
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.