



Recombinant Human Microsomal glutathione S-transferase 1 (MGST1)

Product Code	CSB-CF013791HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P10620
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	VDLTQVMDDEVFMAFASYATIILSKMMLMSTATAFYRLTRKVFANPEDCVAFG KGENAKK YLRTDDRVERVRAHLNDLENIIPFLGIGLLYSLSGPD PSTAILHFRLFVGARIYH TIAY LTPLPQPNRALSFFVGYGVTL SMAYRLLKSKLYL
Source	in vitro E.coli expression system
Target Names	MGST1
Protein Names	Recommended name: Microsomal glutathione S-transferase 1 Short name= Microsomal GST-1 EC= 2.5.1.18 Alternative name(s): Microsomal GST-I
Expression Region	2-155
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.