



Recombinant Sheep V-type proton ATPase 16 kDa proteolipid subunit (ATP6V0C)

Product Code	CSB-CF002389SH
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O18882
Product Type	Transmembrane Protein
Immunogen Species	Ovis aries (Sheep)
Sequence	MSEAKSGPEYASFFAVMGASAAMVFSALGAAYGTAKSGTGIAAMSVMRPEMI MKSIIIPVV MAGIIAIYGLVVAVLIANSLNDGISLYRSFLQLGAGLSVGLSGLAAGFAIGIVGDA GVRG TAQQPRLFVGMILILIFAEVLGLYGLIVALILSTK
Source	in vitro E.coli expression system
Target Names	ATP6V0C
Protein Names	Recommended name: V-type proton ATPase 16 kDa proteolipid subunit Short name= V-ATPase 16 kDa proteolipid subunit Alternative name(s): Vacuolar proton pump 16 kDa proteolipid subunit
Expression Region	1-155
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.