



Recombinant Escherichia coli O157:H7 UPF0442 protein yjjB (yjjB)

Product Code	CSB-CF477944EOE
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B5Z4P6
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli O157:H7 (strain EC4115 / EHEC)
Sequence	MGVIEFLLALAQDMILAAIPAVGFAMVFNVPVRLRWCALLGSIGHGSRMILMT SGLNIE WSTFMASMLGGTIGIQWSRWYLAHPKVFTVA AVIPMFPGISAYTAMISAVKISQ LGYSEP LMITLLTNFLTASSIVGALSIGLSIPGLWLRYKRPRV
Source	in vitro E.coli expression system
Target Names	yjjB
Protein Names	Recommended name: UPF0442 protein yjjB
Expression Region	1-157
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.