



Recombinant Human Transmembrane protein 50A (TMEM50A)

Product Code	CSB-CF023849HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O95807
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MSGFLEGLRCSECIDWGEKRNTIASIAAGVLFFFTGWWIIIDAAVIYPTMKDFNHS YHACG VIATIAFLMINAVSNGQVRGDSYSEGCLGQTGARIWLFVGFMLAFGSLIASMWI LFGGYV AKEKDIVYPGIAVFFQNAFIFFGGLVFKFGRTEDLWQ
Source	in vitro E.coli expression system
Target Names	TMEM50A
Protein Names	Recommended name: Transmembrane protein 50A Alternative name(s): Small membrane protein 1
Expression Region	1-157
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.