



Recombinant Human Transmembrane protein 42 (TMEM42)

Product Code	CSB-CF023839HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q69YG0
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MAERP GPPGGAVSATAYPDTPAEFP PHLQAGAMRRRFWGVFNCLCAGAFGA LAAASAKLA FGSEVSMGLCVLGIIVMASTNSLMWTFFSRGLSFSMSSAIASVTVTFSNILLSSA FLGYVL YGECQEVLLWWGGVFLILCGLTLIHRKLPPTWKPLPHKQQ
Source	in vitro E.coli expression system
Target Names	TMEM42
Protein Names	Recommended name: Transmembrane protein 42
Expression Region	1-159
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.