



Recombinant Megalops atlanticus Cytochrome c oxidase subunit 1 (mt-co1)

Product Code	CSB-CF015072MDE
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P29648
Product Type	Transmembrane Protein
Immunogen Species	Megalops atlanticus (Tarpon) (Clupea gigantea)
Sequence	YQHLEWFFGHPEVYILILPGFGMISHIVAYYAGKKEPFGYMGMVWAMMAIGLL GFIVWAH HMFTVGMVDVTRAYFTSATMIIAIP TGVKVFSWLATLHGGSIKWDTPLLWALGF IFLFTV GGLTGIVLANSSIDIVLHDTYYVVAHFHYVLSMGAVFAIM
Source	in vitro E.coli expression system
Target Names	mt-co1
Protein Names	Recommended name: Cytochrome c oxidase subunit 1 EC= 1.9.3.1 Alternative name(s): Cytochrome c oxidase polypeptide I
Expression Region	1-160
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.