



Recombinant *Macaca fascicularis* Arachidonate 5-lipoxygenase-activating protein (ALOX5AP)

Product Code	CSB-CF001625MOV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q2PG08
Product Type	Transmembrane Protein
Immunogen Species	<i>Macaca fascicularis</i> (Crab-eating macaque) (<i>Cynomolgus</i> monkey)
Sequence	MDQETVGNVLLAIIVTLISVVQNGFFAHKVEHESRTQNGRSFQRTGTLAFERV YTANQNC VDAYPTFLAVLWSAGLLCSQVPAAFAGLMYLLVRQKYFVGYLGERTQSTPGYI FGKRIIL FLFLMSVAGIFNYLIFFFGSDFENYIKTVTTTISPLLLIP
Source	in vitro <i>E.coli</i> expression system
Target Names	ALOX5AP
Protein Names	Recommended name: Arachidonate 5-lipoxygenase-activating protein
Expression Region	1-161
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.