



Recombinant Heme transporter hrg-6 (hrG-6)

Product Code	CSB-CF771998CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q7YTM7
Product Type	Transmembrane Protein
Immunogen Species	Caenorhabditis elegans
Sequence	MRWSLCCHINVRIAYTICGIIIGLFWACVYIFAWKNWVALAACLTATAFAFETFY FYLSV KRDTILNWKSQTFQVLFWINLIVGFLSIGGMITAIVLAATKHQGVSNKDQHGLN WWSTAT WFLVMLKWTWQNAFIARYYGKLLTKSIIHPEEPDDPSTWKF
Source	in vitro E.coli expression system
Target Names	hrG-6
Protein Names	Recommended name: Heme transporter hrg-6 Alternative name(s): Heme-responsive gene 6 protein Short name= CeHRG-6
Expression Region	1-161
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.