



Recombinant *Oryza sativa* subsp. *indica* CASP-like protein BLE3 (BLE3)

Product Code	CSB-CF387251OFF
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A2Y2B7
Product Type	Transmembrane Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>indica</i> (Rice)
Sequence	MAKVHRLMNAVLRRLAAAAAATAAVVMVTSRETTTSFFGIQMEAKYSYTPSFIFF VVAYAV AAAYSLLVLAVPAGSALSRLALTTDVVLGMVLAGAVASAGAISDIKNGNSHAG WLPVCG QIHAYCNHVMAALIAGFVALAVHFVVVMYSLHIVTDVICPCH
Source	in vitro E.coli expression system
Target Names	BLE3
Protein Names	Recommended name: CASP-like protein BLE3 Alternative name(s): Protein brassinolide-enhanced BLE3 Short name= OsBLE3 Short name= Protein BL-enhanced 3
Expression Region	1-162
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.