



Recombinant Schizosaccharomyces pombe Mitochondrial import inner membrane translocase subunit tim17 (tim17)

Product Code	CSB-CF310827SXV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P87130
Product Type	Transmembrane Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Sequence	MASADHTRDPCPYVILNDFGAAFSMGTIGGAIWHSIKGWRNSPPGKRISAIAA AKTRAP VLGGNFGVWGGFLSTFDCAVKGVRRKEDPWNIIAGFFTGGALAVRGGWRA TRNGAIGCA CILAVFEGLGIALGRMNAEYNRPVAPVIPDAPASGSTSAAPAAV
Source	in vitro E.coli expression system
Target Names	tim17
Protein Names	Recommended name: Mitochondrial import inner membrane translocase subunit tim17 Alternative name(s): Mitochondrial protein import protein 2
Expression Region	1-164
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.