



Recombinant *Sulfolobus islandicus* UPF0290 protein LS215_1460 (LS215_1460)

Product Code	CSB-CF502509FPB
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	C3MQ04
Product Type	Transmembrane Protein
Immunogen Species	<i>Sulfolobus islandicus</i> (strain L.S.2.15 / Lassen #1)
Sequence	MSIAYDLLLLSILYLPFAVANGSGPFIKRGTPIDFGKNFVDGRRLFGDGKTFEGLI VALT FGTTVGVIIISKFFTAEWTLISFLESFLFAMIGDMIGAFIKRRLGIPRGGRVGLDQL DFVL GASLILVLMRVNITWYQFLFICGLAFFLHQGTNYVAYLLKIKNVPW
Source	in vitro E.coli expression system
Target Names	carS
Protein Names	Recommended name: UPF0290 protein LS215_1460
Expression Region	1-166
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.