



# Recombinant Methanocaldococcus jannaschii Uncharacterized protein MJ1249.1 (MJ1249.1)

<b>Product Code</b>	CSB-CF305803MRU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P81230
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Methanocaldococcus jannaschii (strain ATCC 43067 / DSM 2661 / JAL-1 / JCM 10045 / NBRC 100440) (Methanococcus jannaschii)
<b>Sequence</b>	MARGINPVLFKISFLLIILVSLILSLFYYNFLFAFLLSFILFGITWAYCYIKIESTDKGK L LYEIKRPGIETLRFLFILMIISVFIKSLLSNSFFPYISFLLSNLILGLVLFDDYILGNP TIKFYEKGVVFDRAVAFYNWEELDIKEDEGYLKIKIKIYYPKEIMYKK
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	MJ1249.1
<b>Protein Names</b>	Recommended name: Uncharacterized protein MJ1249.1
<b>Expression Region</b>	1-166
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.