



Recombinant *Oryza sativa* subsp. *indica* CASP-like protein Osl_30044 (Osl_30044)

Product Code	CSB-CF383748OFF
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A2YXI1
Product Type	Transmembrane Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>indica</i> (Rice)
Sequence	MVELESQEAVTVASTADIAVDVSLRLLAAATSLASAVVVAANHQQRWGVRVDF TLFQVWI GFVAVNLVCTVYAAATAAAARKAMGRWWLHHADAVVVNLEAAATAGAGAIGS IAMWGNEA SGWYAVCRLYRRYCNAGAAALALSAAVLLLGVACARSRYPKMPPTT
Source	in vitro E.coli expression system
Target Names	Osl_30044
Protein Names	Recommended name: CASP-like protein Osl_30044
Expression Region	1-167
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.