



Recombinant *Saccharomyces cerevisiae* Golgi apparatus membrane protein TVP18 (TVP18)

Product Code	CSB-CF312029SVG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q04767
Product Type	Transmembrane Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Sequence	MALSLGQFINVGGMVKDLKSFNFSVYGRWFGYINIILCIALGIANLFHVSGVIAF GIISI IQGLVILFIEIPFLKICPLSDNFIEFIKRFETNGWRCLFYLAMAIQYISIAVMATSLI VVAVGLTISSISYAVAYTKHQEFQNTNIIKNPTDDDFPHEAVVREML
Source	in vitro E.coli expression system
Target Names	TVP18
Protein Names	Recommended name: Golgi apparatus membrane protein TVP18 Alternative name(s): TLG2 compartment vesicle protein of 18 kDa
Expression Region	1-167
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.