



Recombinant *Oryza sativa* subsp. *japonica* CASP-like protein Os08g0536200 (Os08g0536200, LOC_Os08g42430)

Product Code	CSB-CF757910OFG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q6Z1G3
Product Type	Transmembrane Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Sequence	MVELESQEAVTVASTADIAVDVSLRLLAAATSLAAAVVVAANHQQRWGIRVDF TLFQVWI GFVAVNLVCTVYAAATAAAAARKAMGRWWLHHADAVVVNLEAAATAGAGAIG SIAMWGNE ASGWYAVCRLYRRYCNAGAAALALSLAAVLLLGVACARSRYPKMPPTT
Source	in vitro <i>E.coli</i> expression system
Target Names	Os08g0536200
Protein Names	Recommended name: CASP-like protein Os08g0536200
Expression Region	1-168
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.