



Recombinant Populus trichocarpa CASP-like protein POPTRDRAFT_788163 (POPTRDRAFT_788163)

Product Code	CSB-CF497501EXS
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B9NBE5
Product Type	Transmembrane Protein
Immunogen Species	Populus trichocarpa (Western balsam poplar) (Populus balsamifera subsp. trichocarpa)
Sequence	MAAPPAPSMVSRMTALFLRVLTFAFLMVSLVIMTTNTGTIEIGIDEFKVRSKDFY SYRYM LAAIAFGLTYTILQIALTLNHISKRNGAQTSGDGNLVDFYGDKVVS YILATGAAA AFGA TKELKTQLAGLGGDKFFNKGYASASLLLLLGFVCTAILS VFSSYALPKKV
Source	in vitro E.coli expression system
Target Names	POPTR_0010s21230g
Protein Names	Recommended name: CASP-like protein POPTRDRAFT_788163
Expression Region	1-169
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.