



Recombinant Rat Translocator protein (Tspo)

Product Code	CSB-CF025168RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P16257
Product Type	Transmembrane Protein
Immunogen Species	Rattus norvegicus (Rat)
Sequence	MSQSWVPAVGLTLVPSLGGFMGAYFVRGEGLRWYASLQKPSWHPPRWTLA PIWGTLYSAM GYGSYIIWKELGGFTEEAMVPLGLYTGQLALNWAWPPIFFGARQMGWALVDL MLVSGVAT ATTLAWHRVSPPAARLLYPYLAWLAFATMLNYYVWRDNSGRRRGG SRLTE
Source	in vitro E.coli expression system
Target Names	Tspo
Protein Names	Recommended name: Translocator protein Alternative name(s): Mitochondrial benzodiazepine receptor PKBS Peripheral-type benzodiazepine receptor Short name= PBR
Expression Region	1-169
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.