



Recombinant Bovine Translocator protein (TSPO)

Product Code	CSB-CF025168BO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P30535
Product Type	Transmembrane Protein
Immunogen Species	Bos taurus (Bovine)
Sequence	MAPPWVPAVGFTLLPSLGGFLGAQYTRGEGFRWYASLQKPPWHPPRWILAPI WGTLYSAM GYGSYMIWKELGGFSKEAVVPLGLYAGQLALNWAWPPLFFGTRQMGWALVD LLLTGGMAA ATAMAWHQVSPPAACLLYPYLAWLAFAGMLNYRMWQDNQVRRSGRRLSE
Source	in vitro E.coli expression system
Target Names	TSPO
Protein Names	Recommended name: Translocator protein Alternative name(s): Isoquinoline-binding protein Short name= IBP PKBS Peripheral-type benzodiazepine receptor Short name= PBR
Expression Region	1-169
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.