



Recombinant Haemophilus influenzae Uncharacterized transporter HI_0976.1 (HI_0976.1)

Product Code	CSB-CF527075HTA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O86230
Product Type	Transmembrane Protein
Immunogen Species	Haemophilus influenzae (strain ATCC 51907 / DSM 11121 / KW20 / Rd)
Sequence	MAFIGVAILINGGKNNEGIDNISLFGCLLVLSAGIIFA AVL RWTQRVVAKVSTQAY TSVS IVLGTITTL PFTLLL TENWQISLNSTGIAGLLYLAIGCSWLAYWLWNKGLNSVDA NISGV LVALEPLFGILFAVSLLETLSFSAALGITIIMLATLGSTLLPKLLKKS
Source	in vitro E.coli expression system
Target Names	HI_0976.1
Protein Names	Recommended name: Uncharacterized transporter HI_0976.1
Expression Region	1-170
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.