



Recombinant *Oryza sativa* subsp. *japonica* CASP-like protein Os11g0549625 (Os11g0549625, LOC_Os11g34730)

Product Code	CSB-CF646932OFG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q2R2T4
Product Type	Transmembrane Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Sequence	MASSRRTALPSAVLALRLLTLALLAASLAVIAADKLTLDFFGGGLPPKKITFKDVY AYRYV LAIAVIGCAYTLLQIPFVAVSIAKRKRMIGGSENVALLIFADVIFALLVATGAGAG FGL TYDAKSAFGGSKLPGEVVRFFNMAYAAAGLMLLAAAAMALIIMLSIYSLVR
Source	in vitro <i>E.coli</i> expression system
Target Names	Os11g0549625
Protein Names	Recommended name: CASP-like protein Os11g0549625
Expression Region	1-171
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.