



Recombinant *Calidris maritima* NADH-ubiquinone oxidoreductase chain 6 (MT-ND6)

Product Code	CSB-CF015082CAO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P43196
Product Type	Transmembrane Protein
Immunogen Species	<i>Calidris maritima</i> (Purple sandpiper) (<i>Erolia maritima</i>)
Sequence	MTYLVLLGLCFVLGGLAVASNPSPYYGVVGLVLASVAGCGWLLSLGVSFVSL VLFMVYL GGMLVVFVYSVSLAADPFPEAWGDWRVVGYGASFILVVIAGMIVGGLIECWKV GVVTVDS GGMFSVRLDFSGVAMFYSGVGMFLVAGWGLLLTLFVVLELVRGLSRGAIRA V
Source	in vitro E.coli expression system
Target Names	MT-ND6
Protein Names	Recommended name: NADH-ubiquinone oxidoreductase chain 6 EC= 1.6.5.3 Alternative name(s): NADH dehydrogenase subunit 6
Expression Region	1-173
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.