



Recombinant *Cepphus columba* NADH-ubiquinone oxidoreductase chain 6 (MT-ND6)

| | |
|--------------------------|--|
| Product Code | CSB-CF015082CGG |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | P43197 |
| Product Type | Transmembrane Protein |
| Immunogen Species | <i>Cepphus columba</i> (Pigeon guillemot) |
| Sequence | MTYFVVFLGLCFVLGGLAVASNPSPYYGVVGLVLASVVGCGWLLSLGVSFVSL VLFMVYL GGMLVVFVYSVSLAADPFPEAWGDWGVVGYGASFVVVLMVGGVIGGFVGG WDFGVVTVDS VGVFSARLDFSGVAMFYSCGVMFLVAGWGLVLTFLVVLELVRGLSRGAIRAV |
| Source | in vitro E.coli expression system |
| Target Names | MT-ND6 |
| Protein Names | Recommended name: NADH-ubiquinone oxidoreductase chain 6 EC= 1.6.5.3 Alternative name(s): NADH dehydrogenase subunit 6 |
| Expression Region | 1-173 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |