



Recombinant *Saccharomyces cerevisiae* Non-classical export protein 2 homolog 1 (FHN1)

Product Code	CSB-CF346068SVG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P53279
Product Type	Transmembrane Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Sequence	MLSAADNLVRIINAVFLIISIGLISGLIGTQTKHSSRVNFCMFAAVYGLVTDSL YGFAN FWTSLTYPAILLVLDLFLNFIFTFVAATALAVGIRCHSCKNKTYLEQNKIIQGS SSRCHQS QAAVAFFYFSCFLFLIKVTVATMGMMQNGGFGSNTGFSRRRARRQMGIP TISQV
Source	in vitro E.coli expression system
Target Names	FHN1
Protein Names	Recommended name: Non-classical export protein 2 homolog 1 Alternative name(s): Functional homolog of NCE102
Expression Region	1-174
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.