



Recombinant *Bos mutus grunniens* NADH-ubiquinone oxidoreductase chain 6 (MT-ND6)

Product Code	CSB-CF713188BPD
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q5ULB0
Product Type	Transmembrane Protein
Immunogen Species	<i>Bos mutus grunniens</i> (Wild yak) (<i>Bos grunniens</i>)
Sequence	MMLYIVFILSVIFVIGFVGFSSKPSPIYGGLGLIVSGGVGCGIVLNFGGSFLGLMV FLIY LGGMMVVFVGYTTAMATEQYPEIWLLNKAVLGAFITALLMEFFMVYYVLKDKEV EIVFEFN GLGDWVIYDTGDSGFFSEEAMGIAALYSYGTWLVIVTGWSLFIGVVMIMEITRG N
Source	in vitro <i>E.coli</i> expression system
Target Names	MT-ND6
Protein Names	Recommended name: NADH-ubiquinone oxidoreductase chain 6 EC= 1.6.5.3 Alternative name(s): NADH dehydrogenase subunit 6
Expression Region	1-175
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.