



Recombinant Arabidopsis thaliana CASP-like protein At4g25830 (At4g25830)

Product Code	CSB-CF810548DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8L8U9
Product Type	Transmembrane Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Sequence	MVKLRETEVILRLCIVFFLLLTSCLIGLDSQTKEIAYIHKNVSRFRYLLALEAELYIDV VV AAYNLVQLGLGWYNVEQKTSNPKWFSYLLDQTAAYVVFAGTSAAAQHSLLVV TGSRELQW MKWCYKFRTRFCFQMGSAILNYIAAALMVLLSSISAFNLFRLYSPKRFFRFKSSS
Source	in vitro E.coli expression system
Target Names	At4g25830
Protein Names	Recommended name: CASP-like protein At4g25830
Expression Region	1-175
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.