



# Recombinant Mormopterus kalinowskii Cytochrome b (MT-CYB)

<b>Product Code</b>	CSB-CF654988MBAX
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	Q36291
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Mormopterus kalinowskii (Kalinowski's mastiff bat)
<b>Sequence</b>	MTNIRKSHPLIKIVNDAFIDLPA P S N I S S W W N F G S L L G V C L T V Q I L T G L F L A M H Y T S D T A T A F N S V T H I C R D V N Y G W L L R Y L H A N G A S M F F I C L Y L H V G R G L Y Y G S Y T Y T E T W N V G V I L L F A V M A T A F M G Y V L P W G Q M S F W G A T V I T N L L S A I P Y I G T D L V E W I W G G F S V D K A T L T
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	MT-CYB
<b>Protein Names</b>	Recommended name: Cytochrome b Alternative name(s): Complex III subunit 3 Complex III subunit III Cytochrome b-c1 complex subunit 3 Ubiquinol- cytochrome-c reductase complex cytochrome b subunit
<b>Expression Region</b>	1-176
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.