



Recombinant Disulfide bond formation protein B (dsbB)

Product Code	CSB-CF690490SZB
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q54155
Product Type	Transmembrane Protein
Immunogen Species	Shigella flexneri
Sequence	MLRFLNQCSQGRGAWLLMAFTALALELTALWFQHVMLLKPCVLCIYERCALFG VLGAALI GAIAPKTPRLRYVAMVIWLYSAFRGVQLTYEHTMLQLYPSPFATCDFMVRFP EWLPLDKWV PQV FVASGDCAERQWDFLGMEMPQWLLGIFIAYLIVAVLVVISQPFKAKKRD LFRGR
Source	in vitro E.coli expression system
Target Names	dsbB
Protein Names	Recommended name: Disulfide bond formation protein B Alternative name(s): Disulfide oxidoreductase
Expression Region	1-176
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.