



Recombinant Human Vitamin K epoxide reductase complex subunit 1-like protein 1 (VKORC1L1)

Product Code	CSB-CF818673HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8N0U8
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MAAPVLLRVSVPRWERYAVCAAGILLSIYAYHVEREKERDPEHRALCDLG PWVKCSA ALASRWGRGFLLGSIFGKDGVLNQPNSVFGILIFYILQLLLGMTASAVAALILMT SSIMS VVGSLYLAYILYFVLKEFCIICIVTYVLNLLLLIINYKRLVYLNEAWKRQLQPKQD
Source	in vitro E.coli expression system
Target Names	VKORC1L1
Protein Names	Recommended name: Vitamin K epoxide reductase complex subunit 1-like protein 1 Short name= VKORC1-like protein 1
Expression Region	1-176
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.